

Dental Implant Guide System Kit

Patent Number: ZL 2018 2 0803412.2



IG05

- IG05, it is suitable for the implant of single teeth and multiple teeth, determine the location more accurately and quickly.
- Zirconia material, it can be autoclavable at 134°C, with high hardness and high strength, good at wear resistance and oxidation resistance.
- The thermoplastic film can be softened after hot water immersion to cover and make retention on the adjacent teeth. A better operation position can be positioned by the system. Given the location and quantity of implant, the system can be for single tooth and multiple teeth.

Item: IG05


Specification: 150X100.7X45mm

Packing: 1 set/box, 10 boxes/ctn

- Dental implant surgery, due to the rich jaw nerve and blood vessels and individual differences among different patients, it is so hard for dentists to control the position and angles during their drilling in the jaw. And peripheral tissue can be unnecessary damaged. Including: 2 pcs zirconium drills, guide rod 14 pcs, 3 pcs of inner circle, 3 pcs outer circle, partial implant guide 4 pcs, implant guide 12 pcs.



Zirconium drill



Guide rod



Circle



Partial implant guide




Implant guide


Instructions for single teeth implant guide plate:



- Cut the plaster model in central position where needs to do implanting.



- Draw bone topography.




- According to bone topography, draw position and direction of guide rod.



- Take a zirconia drill in kit, grind 1.1mm groove on both side of plaster model, insert the guide rod, and then glue the two half plaster model together.




- Take single teeth implant guide plate, soak and soften it into hot water, after the color becomes transparent, take out.



- Place the single teeth guide plate on the plaster model, and make shape.



- Take down when the plate becomes harder.



- Black inner ring suits for drill below 2.3mm.



- White inner ring suits for drill below 3.0mm.

DISPOSABLE DENTAL AIR WATER SYRINGE TIPS



Disposable Dental Air Water Syringe Tips



AWT01


■ AWT01,universal three-in-one Air Water Syringe Tips. The inner tube is made of stainless steel and the outer tube is made of white plastic.

Item: AWT01
Specification: 84X ϕ 3.87mm
Packing: 250 pcs/box,30 boxes/ctn
Carton Size: 290X260X230mm
Weight: 11.6 KGS



White

Disposable Dental Air Water Syringe Tips



AWT02


■ AWT02,white plastic material. Independent plastic tube structure without core. There are two modes with Syringe to use: water spray and water vapor.

Item: AWT02
Specification: 80X ϕ 3.87mm
Packing: 250 pcs/box,30 boxes/ctn
Carton Size: 290X260X230mm
Weight: 7 KGS



White

Disposable Dental Air Water Syringe Tips



AWT03

■ AWT03,universal three-in-one Air Water Syringe Tips. Hard plastic material. The inner tube is colored and the outer tube is transparent plastic

Item: AWT03
Specification: 92X ϕ 3.87mm
Packing: 250 pcs/box,40boxes/ctn
Carton Size: 420X270X230mm
Weight: 13.1 KGS

Light Yellow
Light Orange
Light Green
Light Blue
Light Purple
White




Disposable Dental Air Water Syringe Tips



AWT05

Item: AWT05
Specification: 76X ϕ 3.87mm
Packing: 250 pcs/box,40 boxes/ctn
Carton Size: 340X270X230mm
Weight: 10 KGS



DAM KIT

