ENDO SYSTEM







Fi-P



0.2 second heating time


High-performance conversion system ensures the efficient cutting





High precision temperature control

Constant temperature prevents scalding





4 hours of endurance

Large capacity battery & special heating system enable 1500 operations with a full charge.





0.35mm pre-bend tip

Thinner tips can penetrate deep into the narrow root canal, providing you a better treatment effect.






4 working temperatures

150°C、180°C、200°C、230°C can match gutta-percha of different brands on the market.





WP5508L





Quick heating within 15 seconds

reduces time for waiting and increase your working efficiency.





0.5mm injection needle

Easy to reach the apical 3mm of canal brings better filling.





WP3504 WP4004

WP4504

WP5004

WP5506

360 degree rotatable injection needle

0.50mm

0.55mm

reduces time for waiting and increase your working efficiency.

0.04

0.06



Dual screen display

Unique double displays design enables convenient operations in different angles and teeth positions.



Taper	Size	Length	Taper	Size	Length
20G 22mm	20G	22mm	23G 24mm	23G	24mm
20G 24mm	20G	24mm	23G 28mm	23G	28mm
20G 28mm	20G	28mm	25G 24mm	25G	24mm