

General Consumables

Robotic Pipette Tips

A liquid handling robot is used to automate workflows in life science laboratories. It takes the manual labor out of repeated pipetting with improved accuracy and repeatability. Pipette tips can be made conductive to prevent carryover in the automated pipetting process

Pipette tips available are compatible with:

- Tecan
- Hamiliton
- OpenTron
- · Agilent Bravo
- Beckman
- PerkinElmer
- Olympus







OpenTron OT-2 Compatible Tips

| Cat No | $Max\;Volume\;(\mu L)$ | Specifications | Packagaing |
|---------------|------------------------|----------------------------------|------------------------|
| TF20-RB-CLS | 20 | Low retention, filtered, sterile | 96pc/rack, 50racks/ctn |
| TF200-RB-CLS | 200 | Low retention, filtered, sterile | 96pc/rack, 50racks/ctn |
| TF1000-RB-CLS | 1000 | Low retention, filtered, sterile | 96pc/rack, 50racks/ctn |

Agilent Bravo Compatible Tips

| Cat No | Max Volume (μL) | Specifications | Packaging |
|--------------|-----------------|---------------------------|----------------------------------|
| MF-RE-384 | 30 | Non-filtered, Non-sterile | 384pc/rack, 10rack/box, 5box/ctn |
| MF-RE-384 | 30 | Filtered, Non-sterile | 384pc/rack, 10rack/box, 5box/ctn |
| MF-RE-384 | 70 | Non-filtered, Non-sterile | 384pc/rack, 10rack/box, 5box/ctn |
| MF-RE-384 | 70 | Filtered, Non-sterile | 384pc/rack, 10rack/box, 5box/ctn |
| RG-250-R-CS | 250 | Non-filtered, Sterile | 384pc/rack, 10rack/box, 5box/ctn |
| RGF-250-R-CS | 250 | Filtered, Sterile | 384pc/rack, 10rack/box, 5box/ctn |



General Consumables

Tecan Compatible Tips

| Land State . | Cat No | Volume (µL) | Specifications | Packaging |
|--------------|----------------------|-------------|-----------------------------------|-----------------------|
| | MF-TE-96 | 50 | Conductive, Filtered, Sterile | 96pc/rack, 24rack/ctn |
| | MF-TE-96 | 200 | Conductive, Filtered, Sterile | 96pc/rack, 24rack/ctn |
| | MF-TE-96 | 1000 | Conductive, Filtered, Sterile | 96pc/rack, 24rack/ctn |
| - diffine - | MF-TE-96 (Wide Bore) | 1000 | Conductive, Non-filtered, Sterile | 96pc/rack, 24rack/ctn |
| | MF-TE-96 (Wide Bore) | 1000 | Conductive, Filtered, Sterile | 96pc/rack, 24rack/ctn |
| | RTF-50-H-ES | 50 | Conductive, Filtered, Sterile | 96pc/rack, 36rack/ctn |
| | RTF-200-H-ES | 200 | Conductive, Filtered, Sterile | 96pc/rack, 36rack/ctn |
| | RTF-1000-H-ES | 1000 | Conductive, Filtered, Sterile | 96pc/rack, 36rack/ctn |
| | RTF-50-H-CS | 50 | Filtered, Sterile | 96pc/rack, 36rack/ctn |
| | RTF-200-H-CS | 200 | Filtered, Sterile | 96pc/rack, 36rack/ctn |
| | RTF-1000-H-CS | 1000 | Filtered Sterile | 96nc/rack 36rack/ctn |

Wide Bore Tecan Tips

Beckman Compatible Tips

| Cat No | Max Volume (µL) | Specifications | Packaging |
|---------------|-----------------|-----------------------|---------------------------------|
| RB-50-R-CS | 50 | Non-filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RBF-50-R-CS | 50 | Filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RB-250-R-CS | 250 | Non-filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RBF-50-R-CS | 250 | Filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RB-1000-R-CS | 1000 | Non-filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RBF-1000-R-CS | 1000 | Filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |

Hamilton Compatible Tips

| Cat No | Max Volume (μL) | Specifications | Packaging |
|--------------|-----------------|-----------------------------------|-----------------------|
| RH50-H-ES | 50 | Non-filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RHF50-H-ES | 50 | Filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RH50-H-ES | 300 | Non-filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RHF50-H-ES | 300 | Filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RH1000-H-ES | 1000 | Non-filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RHF1000-H-ES | 1000 | Filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RH50-H-CS | 300 | Non-filtered, Sterile | 96pc/rack, 36rack/ctn |
| RHF50-H-CS | 300 | Non-filtered, Sterile | 96pc/rack, 36rack/ctn |
| RH1000-H-CS | 1000 | Non-filtered, Sterile | 96pc/rack, 36rack/ctn |
| RHF1000-H-CS | 1000 | Non-filtered, Sterile | 96pc/rack, 36rack/ctn |
| | | | |