



Item: SG04
Packing: 12 pcs/box,18 boxes/ctn
Carton Size: 500X360X430mm
Weight: 10 KGS

Yellow ☒
Clear ☐



Safety Glasses



Item: SG06
Packing: 12 pcs/box,18 boxes/ctn
Carton Size: 500X360X430mm
Weight: 10 KGS



Disposable Eyeshield(Anti-fog)



Item: DE01
Packing: 1 frame+1 shield/set,5 sets/bag,100 bags/ctn
Carton Size: 430X270X240mm
Weight: 5.4 KGS
Remarks: (shields)10 pcs/bag,300 bags/ctn

Yellow ☒
Green ☐
Blue ☐
Purple ☐



Protective Light Hand Shield



Item: YK03
Specification: 127X240X3mm
Packing: 1 pc/bag,100 bag/ctn
Carton Size: 290X260X230mm
Weight: 8.3 KGS

SECTIONAL CONTOURED MATRICES KIT



Dental Sectional Contoured Matrices Kit



SC1001

Item: SC1001/SC1002
Name: Dental Sectional Contoured Matrices Kit
Packing: (SC1001)1set/box,10 boxes/ctn
(SC1002)1set/box,40 boxes/ctn



SC1002

SC1001 Contains:

 PW01 20pcs PW02 20 pcs PW03 20 pcs Dental Wedge	 SC1201L 50 pcs Ordinary Matrix Band	 SC1202M 50 pcs Ordinary Matrix Band	 SC1203S 50 pcs Ordinary Matrix Band	 SC0601L 10 pcs SC0601S 10 pcs Resin Matrix Band SC0301 6 pcs Silicone Clearance Wedge	 SC08 3 pcs Adhesive Tip Applicator
 SC0206 1pc SC0205 1pc SC0204 1pc	 SC1204L 20 pcs Subgingival Matrix Band	 SC1205M 20 pcs Subgingival Matrix Band	 SC1206S 20 pcs Subgingival Matrix Band	 WedgeGuard WG01 15 pcs Dental Wedge with Matrix PW05(left)15 pcs (right)15 pcs Left Right Dental Wedge with Matrix PW06 10 pcs	 SC0702 1pc Forceps

SC1002 Contains:

 SC1201L 15 pcs Ordinary Matrix Band	 SC1204L 5 pcs Subgingival Matrix Band	 SC0206 SC0205 SC0204 Clamping Ring	 SC1203S 15 pcs Ordinary Matrix Band	 SC1206S 5 pcs Subgingival Matrix Band
 PW01 4pcs PW02 4pcs PW03 6pcs Dental Wedge	 SC1202M 15 pcs Ordinary Matrix Band	 SC1205M 5 pcs Subgingival Matrix Band	 WG01 3pcs WedgeGuard	 SC0301 6 pcs Silicone Clearance Wedge

Dental Sectional Contoured Matrices Kit



SC01

Item: SC01/SC0104
Name: Dental Sectional Contoured Matrices Kit
Packing: (SC01)1set/box,10 boxes/ctn
(SC0104)1set/box,40 boxes/ctn
Carton Size: (SC01)310X250X270mm
(SC0104)355X345X265mm
Weight: (SC01)4.9 KGS



SC01 Contains:

 PW01 20pcs PW02 20 pcs PW03 20 pcs Dental Wedge	 SC0501L 50 pcs Ordinary Matrix Band	 SC0502M 50 pcs Ordinary Matrix Band	 SC0503S 50 pcs Ordinary Matrix Band	 SC0601L 10 pcs SC0601S 10 pcs Resin Matrix Band	 SC08 3 pcs Adhesive Tip Applicator
 SC0201 2pcs SC0202 2pcs SC0203 2pcs Clamping Ring	 SC0504L 20 pcs Subgingival Matrix Band	 SC0505M 20 pcs Subgingival Matrix Band	 SC0506S 20 pcs Subgingival Matrix Band	 SC0301 6 pcs Silicone Clearance Wedge	 SC07 1pc Forceps

SC0104 Contains:

 SC0501L 15 pcs Ordinary Matrix Band	 SC0504L 5 pcs Subgingival Matrix Band	 SC0206 SC0205 SC0204 Clamping Ring	 SC0503S 15 pcs Ordinary Matrix Band	 SC0506S 5 pcs Subgingival Matrix Band
 PW01 4pcs PW02 4pcs PW03 6pcs Dental Wedge	 SC0502M 15 pcs Ordinary Matrix Band	 SC0505M 5 pcs Subgingival Matrix Band	 WG01 3pcs WedgeGuard	 SC0301 6 pcs Silicone Clearance Wedge







Dental Sectional Contoured Matrices Kit



SC0103

Item: SC0103
Packing: 1 set/box,100 boxes/ctn
Carton Size: 420X270X230mm
Weight: 7 KGS

The Kit Contains:

Ordinary Matrix Band (L) SC0501L 20 pcs Subgingival Matrix Band (L) SC0504L 10 pcs 	Ordinary Matrix Band (M) SC0502M 25 pcs Subgingival Matrix Band (M) SC0505M 15 pcs 	Ordinary Matrix Band (S) SC0503S 20 pcs Subgingival Matrix Band (S) SC0506S 10 pcs 	Clamping Ring SC0202 1 pc SC0203 1 pc 	Dental Wedge PW01 8 pcs PW02 8 pcs Size: L, M, S, SS 	WedgeGuard WG01 9 pcs Dental Wedge with Matrix PW05 12 pcs, PW06 4 pcs 
--	--	--	---	---	--

Clamping Ring

NEW




It consists of U-shaped alloy spring clips and plastic chucks. The materials used are titanium-nickel memory alloy and PC plastic, which are super-elastic, memory, corrosion-resistant and biocompatible. The clamping force is weak, medium and strong respectively.

Advantage:

- U-shaped memory alloy clamping ring to ensure the strength of the ring;
- Soft and hard plastic clamping heads, play a anti-slip effect;
- Head covered with transparent material, which does not affect light curing restoration;
- The design of the V-shaped clamping head can fix the two ends of the interdental space and make full contact with the wedge for more accurate positioning.

Item: SC0204/SC0205/SC0206
Specification: 25.5X22mm
Packing: 1 pc/box, 30 boxes/ctn
Carton Size: 210X150X140mm

Red (Strong)
Blue (Medium)
Green (Weak)



Clamping Ring



The concave type design of the silicone pads conforms with the anatomy of the tooth adjacent surface, closely fixed on the matrix bands on it, make the adjacency relations back to normal. The silicone pads circle, vertical circle and triangle all of them can be used to fix the bands.

SC0401 Silicone pad
Packing: 50 pcs/box

Item: SC0201/SC0202/SC0203
Packing: 6 pcs/box, 68 boxes/ctn
Carton Size: 260X170X180mm



Silicone Clearance Wedge

Patent Number: ZL 2015 2 0937846.8




"8" appearance, hollow design, Super elasticity, can be easily put into any diastema. It is used to fix the resin matrix bands, suitable for repairing the smaller cavity of anterior teeth.

Item: SC0301
Packing: 20 pcs/box, 68 boxes/ctn
Carton Size: 260X170X180mm



Resin Matrix Band



Made of transparent material, without affecting of the resin curing, easily of operation.

SC0601L / SC0602S

Item: SC0601L/SC0602S
Packing: 50 pcs/box, 68 boxes/ctn
Carton Size: 260X170X180mm



Ordinary Matrix Band & Subgingival Matrix Band

Patent Number: ZL 2019 2 1915792.X

NEW



Ordinary Matrix Band

Subgingival Matrix Band



SC1201L



SC1204L



SC1202M



SC1205M



SC1203S



SC1206S

■ Anti-Stick Type: It's coated with Teflon non-stick material, which can prevent the matrix band from sticking to resin material during dental restoration, and make the filling resin surface smoother and reduce the adhesion of plaque.

■ Size for Anti-Stick Type: large purple, medium orange, small green.



Ordinary Matrix Band & Subgingival Matrix Band



Ordinary Matrix Band

Subgingival Matrix Band



■ Halves bean appearance, smooth surface, convenient operation. Concave type design ordinary matrix band and subgingival matrix bands make it easy to recover to the normal adjacency relations.

Ordinary Matrix Band



SC0501L : 13.7X6.8mm
SC0502M : 13X5.6mm
SC0503S : 11.7X4.8mm

Subgingival Matrix Band



SC0504L : 14X9mm
SC0505M : 12.5X6.7mm
SC0506S : 11.5X5mm



Packing: 50 pcs/box, 68 boxes/ctn
Carton Size: 260X170X180mm

Resin Matrix Band



SC09

■ The product is made of polyester film, which can restore the defect tooth shape well. Secured with poly wedge, the matrix band can fix tooth cervical firmly, reduce overhang to gum, and restore proximal contact point.

Item: SC09

Specification: ■ 25X10X0.075mm 100pcs / 100X10X0.075mm 20pcs
■ 25X10X0.05mm 100pcs / 100X10X0.05mm 20pcs

Packing: 1 set/box, 60 boxes/ctn

Carton Size: 360X170X280mm



Forceps



SC0702

Item: SC0702

Packing: 1 pc/box, 50 boxes/ctn



Forceps



SC07

Item: SC07

Packing: 1 pc/box, 50 boxes/ctn

Carton Size: 280X280X220mm

Weight: 4.6 KGS



Adhesive Tip Applicator

Patent Number: ZL 2015 2 0375955.5



SC08


■ Used to stick small device like ceramic veneer, inlay and Orthodontic bracke, special made glue tip with strong stick force but not remians left over, well-designed, convenience to use to avoid getting slipped from the forceps.

Item: SC08

Packing: 20 pcs/box, 120 boxes/ctn

Carton Size: 440X440X250mm

Weight: 6 KGS




WedgeGuard

Patent Number: ZL 2017 2 1625810.1




■ According to the size of adjacent gap,select the appropriate WedgeGuard and insert with forceps or handle from the buccal or lingual side between the desired preparation teeth and adjacent teeth.Then adjust the position to fit the tooth better.Later,the tooth preparation can be performed.



Item: WG01
Packing: 50 pcs/box,160 boxes/ctn
Carton Size: 280X280X220mm
Weight: 9.6 KGS



Item No	Color	A(±0.1mm)	B(±0.1mm)	C(±0.1mm)	Thickness
WG01L	Light Purple	2.8 mm	1.9 mm	2.4 mm	0.06 mm
WG01M	Light Orange	2.0 mm	1.3 mm	1.7 mm	0.06 mm
WG01S	Light Green	1.9 mm	0.9 mm	1.4 mm	0.06 mm



Dental Wedge with Matrix




■ This product is composed of Wedges and Matrix Bands,which can not be separated.The product is made of transparent plastic,which is convenient for curing and forming by light curing machine after filling with resin.It can effectively shorten the curing time and achieve excellent molding effect.The Dental Wedge with Matrix is designed to touch the patient' s dental papilla in a rounded way.The outside design can effectively prevent slip,and feel good.The operation is convenient and quick.This product is an important auxiliary tool for the dentist in the key process of repair and healing.



Item: PW05/PW06
Packing: 50(L+R) pcs/box,100 boxes/ctn
Carton Size: 290X260X230mm


Item No	Specification	Color
PW05L(L/R)	19.5x6.7mm	Blue, Clear
PW05M(L/R)	17.6x6.8mm	
PW05S(L/R)	15.6x6.8mm	
PW06L	17.7x8.8mm	
PW06S	17.7x8.2mm	



Common Plastic Wedge



■ PW01:The design of the round tip protect gingival from damage,the square tail is very convenient to hold .



Item: PW01
Packing: 100 pcs/bottle,16 bottles/box,36boxes/ctn
Carton Size: 430X430X290mm
Weight: (L,M)12 KGS,(S,SS)10 KGS

Item No	Color	Specification
PW01L	Purple	16X2 . 6X2 . 2mm
PW01M	Yellow	14X2 . 2X2mm
PW01S	Green	12X1 . 9X1 . 9mm
PW01SS	Blue	10X1 . 6X1 . 5mm




Plastic Wedge with Hole

Patent Number: ZL 2014 2 0321698.2



■ PW02 : Plastic wedge,round stern with a hole,easily of being taken,round edge with unattached,make the wedge more flexibility to be full filled in different teeth cracks .



Item: PW02
Packing: 100 pcs/bottle,16 bottles/box,36boxes/ctn
Carton Size: 430X430X290mm
Weight: (L,M)11 KGS,(S,SS)9 KGS

Item No	Color	Specification
PW02L	Purple	17 . 8X2 . 8mm
PW02M	Yellow	17 . 8X2 . 45mm
PW02S	Green	17 . 8X2 . 1mm
PW02SS	Blue	17 . 8X1 . 7mm



Wooden Wedge



■ PW03 : Eco-friendly wooden wedge,round tips protect gingival from damage,smooth round edge, square stern.Six different sizes fits for different requirements .



Item: PW03
Packing: 100 pcs/bottle,16 bottles/box,36boxes/ctn
Carton Size: 430X430X290mm
Weight: 9-11 KGS

Color	L(±0.5mm)	X(±0.1mm)	Y(±0.1mm)	H(±0.5mm)
Orange	11mm	3.5mm	1.5mm	1.3mm
White	14mm	4mm	1.5mm	1.5mm
Green	14mm	2mm	1.5mm	1.5mm
Blue	14mm	4mm	2mm	1.5mm
Yellow	16mm	4mm	1.5mm	1.5mm
Pink	16mm	4.5mm	2mm	1.5mm




Long Handle Wedge

Patent Number: ZL 2015 2 0018098.3



■ PW04 : Long handle plastic wedge,integrated design,just break off the long handle lightly when the wedge is in place .



Item: PW04
Packing: 100 pcs/bag,100 bags/ctn
Carton Size: 290X260X230mm
Weight: (L)6 KGS,(M)5.9 KGS,(S)5.8 KGS

Item No	Color	Specification
PW04L	Orange	81X2 . 7mm
PW04M	Blue	80X2 . 6mm
PW04S	Yellow	75X2 . 1mm



Autoclavable Plastic Tray



■ It can be filled with different kinds of medical instruments, against contamination.

Item: APT01
Specification: 336x240x23mm
Packing: 40 pcs/ctn
Carton Size: 370X270X245mm
Weight: 10 KGS

Green ☒
Blue ☒



Disposable Plastic Tray



Item: DPT01
Specification: 190X148X15mm
Packing: 400 pcs/ctn
Carton Size: 310X210X195mm
Weight: 2.65 KGS

White ☐



Disposable Plastic Tray



Item: DPT02
Specification: 190X148X15mm
Packing: 400 pcs/ctn
Carton Size: 310X210X195mm
Weight: 2.68 KGS

White ☐



Disposable Plastic Tray



Item: DPT03
Specification: 289X190X17mm
Packing: 400 pcs/ctn
Carton Size: 400X305X140mm
Weight: 5.06 KGS

White ☐



Disposable Plastic Tray



Item: DPT04
Specification: 284X186X17mm
Packing: 400 pcs/ctn
Carton Size: 400X305X140mm
Weight: 5.1 KGS

White ☐



Autoclavable Instrument Kit



Item: AK01
Packing: 60 pcs/ctn
Carton Size: 520X410X190mm
Weight: 5.9 KGS-8 KGS

Compule Dispesnser Gun



■ Compule dispenser gun is used for conveying ceramic veneer bonding resin.
■ It can be used repeatedly, need to do heat sterilization at 134°C before using. (Note: Please disassemble the product as below firstly, then do sterilization.)

Item: CG01
Specification: 150X100mm
Packing: 160 pcs/ctn
Carton Size: 420X320X290mm
Weight: 9 KGS



Lift the Clip



Apart Firstly, then Sterile it

Compule Dispesnser Gun Bullets



■ The bullrts material can avoid the light, and is a safe injection resin material. Bullet is disposable, and can prevent cross infection.

Item: CG02
Packing: 100 pcs/bag, 100 bags/ctn
Carton Size: 410X270X230mm
Weight: 5 KGS

Black ☒



Alcohol Lamp



Item: AL01
Capacity: 100ml
Packing: 1 pc/box, 60 boxes/ctn
Carton Size: 520X430X340mm
Weight: 13.9 KGS



Neaten Case



■ For technical use neaten case, the clamping piece on the slice can be used to fix paper, color variety, convenient for distinguishment.

Item: NC01
Specification: 154X138X205mm
Packing: 50 sets/ctn
Carton Size: 430X330X300mm
Weight: 7.7 KGS

Yellow ☒
Purple ☒
Orange ☒
Blue ☒
Red ☒
Black ☒
White ☐



Curved Utility Syringe



Item: CS01
Capacity: 12ml
Packing: 50 pcs/bag, 25 bags/ctn
Carton Size: 500X360X430mm
Weight: 12.4 KGS



Straight Syringe



Item: CS02
Capacity: 5ml
Packing: 50 pcs/bag, 20 bags/ctn
Carton Size: 510X345X370mm
Weight: 8 KGS



Disposable Dental Examination Kit(8in1)



DK04

■ It is used for oral diagnosis and treatment,through strict disinfection and sterilization of ethylene oxide,which can prevent various diseases caused by cross infection.Including:mouth mirror,probe,forceps,bibs,facial tissue,tray,cotton balls,saliva ejector.

Item: DK04
Packing: 200 sets/ctn
Carton Size: 600X390X310mm
Weight: 10 KGS

Disposable Dental Examination Kit(7in1)



DK03

■ Ideal for the diagnosis and treatment of dental,avoid of all kinds of diseases because of cross-infection.Sterilized by ethylene oxide strictly,two years valid date.Including:mouth mirror,probe,forceps,bibs,facial tissue,tray,cotton balls.

Item: DK03
Packing: 200 sets/ctn
Carton Size: 600X390X310mm
Weight: 9 KGS

Disposable Dental Examination Kit(3in1)



DK01A

Item: DK01A
Packing: 1 set/pouch,20 pouches/bag,25 bags/ctn
Carton Size: 420X330X310mm
Weight: 10.5 KGS

X-Ray Film Holder




XH01

■ Angled with both ends working to take both anterior and posterior films. Help finishing the technique photograph of corecting the angular bisector.

Item: XH01
Specification: 126x25x11mm
Packing: 4 pcs/box,200 boxes/ctn
Carton Size: 430X330X310mm
Weight: 16 KGS



Universal Barrier Film



UF

■ Universal barrier film protect the equipment against other pollution,and the smart design make people use it more easily.Total 1200 sheets,length of each sheet:10 cm.

Item: UF
Packing: 1 pc/box,8 boxes/ctn
Carton Size: 360X350X280mm
Weight: 10 KGS
Remark: With Edge/Without Edge

Blue ☒
White ☐



Barrier Film Dispenser



TH06B

Item: TH06B
Specification: 242x137x180mm
Packing: 1 pc/box,12 boxes/ctn
Carton Size: 530X380X480mm
Weight: 9 KGS



Dental Bib Clip



■ Silicone material for the colorized part,stainless material for the head part.

Item: QI15C
Specification: Length of silicone material:397mm;
Head part length:51mm
Packing: 1 pc/pouch,5 pouches/bag,100 bags/ctn
Carton Size: 400X400X370mm
Weight: 10 KGS

Light Yellow ☒
Light Orange ☐
Light Green ☐
Light Blue ☐
Light Purple ☐



Resin Syringe Holder




Item: TH02L/TH02S
Specification: (TH02L)203X108X130mm,
(TH02S)102X108X130mm



Dental Bib Box




Item: TH07
Specification: 175X135X200mm



Multi-function Holder



Item: TH04MFL/TH04MFS
Specification: (TH04MF1)192X115X141mm,
(TH04MF2)112X112X141mm



Multi-function Holder



Item: TH04MF7
Specification: 210x160x95mm

