




General Consumables

Universal Pipette Tips

Manufactured from high-quality virgin polypropylene to fit a wide variety of single and multi-channel pipettes.

- Low Retention available
- Aerosol filter available
- Extended tips available
- Racked and bagged packaging
- DNase- and RNase- free
- Sterilized via gamma irradiation



Low Retention, Filtered, Racked and Sterile Pipette Tips

| Cat No | Max Volume (µL) | Length (mm) | Specifications | Packagaing |
|---------------|-----------------|-------------|----------------------------------|---------------------------------|
| MF-TR-96-10 | 10 | 46 | Low retention, filtered, sterile | 96pc/rack, 10rack/box, 5box/ctn |
| MF-TR-96-20 | 20 | 60 | Low retention, filtered, sterile | 96pc/rack, 10rack/box, 5box/ctn |
| MF-TR-96-100 | 100 | 60 | Low retention, filtered, sterile | 96pc/rack, 10rack/box, 5box/ctn |
| MF-TR-96-200 | 200 | 60 | Low retention, filtered, sterile | 96pc/rack, 10rack/box, 5box/ctn |
| MF-TR-96-300 | 200 | 60 | Low retention, filtered, sterile | 96pc/rack, 10rack/box, 5box/ctn |
| MF-TR-96-1000 | 1000 | 78 | Low retention, filtered, sterile | 96pc/rack, 10rack/box, 5box/ctn |
| MF-TR-96-1250 | 1000 | 102 | Low retention, filtered, sterile | 96pc/rack, 10rack/box, 5box/ctn |


Enquire for bagged, non filtered, non sterile tips or short 33mm 10uL pipette tips



General Consumables

Eppendorf, Satorius Compatible Tips

| Cat No | Max Volume (μL) | Specifications | Packaging |
|--------------|-----------------|----------------------------------|-----------------------|
| TLF10-R-CS | 10 | Low Retention, Filtered, Sterile | 96pc/rack, 50rack/ctn |
| TLF200-R-CS | 200 | Low Retention, Filtered, Sterile | 96pc/rack, 50rack/ctn |
| TLF1000-R-CS | 1000 | Low Retention, Filtered, Sterile | 96pc/rack, 50rack/ctn |
| TL10-B-C | 10 | Low Retention, bagged | 1000pc/bag, 10bag/ctn |
| TL200-B-C | 200 | Low Retention, bagged | 1000pc/bag, 10bag/ctn |
| TL1000-B-C | 1000 | Low Retention, bagged | 1000pc/bag, 10bag/ctn |



39mm, 46mm 10μL Tips available




Eppendorf, Satorius Low Retention Quality

| Cat No | Water | 75% Ethanol | X-100 | SDS |
|--------------|--------|-------------|-------|-------|
| TL-10-R-CS | ±0.1% | ±0.6% | ±0.4% | ±0.6% |
| TL-200-R-CS | ±0.1% | ±0.3% | ±1.1% | ±1.6% |
| TL-1000-R-CS | ±0.02% | ±0.6 | ±1.5 | ±2.1 |



Rainin Compatible Tips

| Cat No | Volume (μL) | Specifications | Packaging |
|----------------|-------------|----------------------------------|-----------------------|
| RAF-10-R-CLS | 10 | Low retention, Filtered, Sterile | 96pc/rack, 50rack/ctn |
| RAF-20-R-CLS | 20 | Low retention, Filtered, Sterile | 96pc/rack, 50rack/ctn |
| RAF-200-R-CLS | 200 | Low retention, Filtered, Sterile | 96pc/rack, 50rack/ctn |
| RAF-250-R-CLS | 250 | Low retention, Filtered, Sterile | 96pc/rack, 50rack/ctn |
| RAF-1000-R-CLS | 1000 | Low retention, Filtered, Sterile | 96pc/rack, 50rack/ctn |





General Consumables

Reload Systems

- Pipette tips in Reload System come pre-loaded on cards for quick and easy loading into empty racks.
- Reload System for pipette tips saves time in refilling empty pipette racks.
- Non-filtered, nonsterile

| Cat No | Volume (μL) | Specifications | Packaging |
|--------------|-------------|----------------|------------------------|
| MF-96RL-10 | 10 | Low retention | 960pc/tray, 10tray/ctn |
| MF-96RL-200 | 200 | Low retention | 960pc/tray, 10tray/ctn |
| MF-96RL-300 | 300 | Low retention | 960pc/tray, 10tray/ctn |
| MF-96RL-1000 | 1000 | Low retention | 960pc/tray, 10tray/ctn |
| MF-96RL-1250 | 1250 | Low retention | 960pc/tray, 10tray/ctn |



Gel Loading Pipette Tip, BS-200-GTB

Gel Loading Pipette Tips are manufactured from high-quality polypropylene material and facilitate precision sample loading into gels. Clear polypropylene tips include long, thin capillary sections. Useful for a wide variety of sample loading applications such as protein separation and other gel loading applications.

- 200uL, bulk packed
- 1000pc/bag, 10bag/ctn
- DNase- and RNase- free
- Nonpyrogenic
- Nonsterile



Reagent Reservoir, BS-PS50,RR

- Easy-to-use loading of multichannel pipettes
- A wide base provides stability and avoids spills
- "V" shaped well design allows maximum reagent recovery
- Polystyrene reservoir are in sterile form, with graduations
- Sterile
- 1pc/bag, 100bag/ctn





General Consumables

Robotic Pipette Tips


A liquid handling robot is used to automate workflows in life science laboratories. It takes the manual labor out of repeated pipetting with improved accuracy and repeatability. Pipette tips can be made conductive to prevent carryover in the automated pipetting process

Pipette tips available are compatible with:

- Tecan
- Hamilton
- OpenTron
- Agilent Bravo
- Beckman
- PerkinElmer
- Olympus



OpenTron OT-2 Compatible Tips



| Cat No | Max Volume (μL) | Specifications | Packaging |
|---------------|-----------------|----------------------------------|------------------------|
| TF20-RB-CLS | 20 | Low retention, filtered, sterile | 96pc/rack, 50racks/ctn |
| TF200-RB-CLS | 200 | Low retention, filtered, sterile | 96pc/rack, 50racks/ctn |
| TF1000-RB-CLS | 1000 | Low retention, filtered, sterile | 96pc/rack, 50racks/ctn |

Agilent Bravo Compatible Tips




| Cat No | Max Volume (μL) | Specifications | Packaging |
|--------------|-----------------|---------------------------|----------------------------------|
| MF-RE-384 | 30 | Non-filtered, Non-sterile | 384pc/rack, 10rack/box, 5box/ctn |
| MF-RE-384 | 30 | Filtered, Non-sterile | 384pc/rack, 10rack/box, 5box/ctn |
| MF-RE-384 | 70 | Non-filtered, Non-sterile | 384pc/rack, 10rack/box, 5box/ctn |
| MF-RE-384 | 70 | Filtered, Non-sterile | 384pc/rack, 10rack/box, 5box/ctn |
| RG-250-R-CS | 250 | Non-filtered, Sterile | 384pc/rack, 10rack/box, 5box/ctn |
| RGF-250-R-CS | 250 | Filtered, Sterile | 384pc/rack, 10rack/box, 5box/ctn |

Enquire for other specifications or packaging



General Consumables

Tecan Compatible Tips



| Cat No | Volume (µL) | Specifications | Packaging |
|----------------------|-------------|-----------------------------------|-----------------------|
| MF-TE-96 | 50 | Conductive, Filtered, Sterile | 96pc/rack, 24rack/ctn |
| MF-TE-96 | 200 | Conductive, Filtered, Sterile | 96pc/rack, 24rack/ctn |
| MF-TE-96 | 1000 | Conductive, Filtered, Sterile | 96pc/rack, 24rack/ctn |
| MF-TE-96 (Wide Bore) | 1000 | Conductive, Non-filtered, Sterile | 96pc/rack, 24rack/ctn |
| MF-TE-96 (Wide Bore) | 1000 | Conductive, Filtered, Sterile | 96pc/rack, 24rack/ctn |
| RTF-50-H-ES | 50 | Conductive, Filtered, Sterile | 96pc/rack, 36rack/ctn |
| RTF-200-H-ES | 200 | Conductive, Filtered, Sterile | 96pc/rack, 36rack/ctn |
| RTF-1000-H-ES | 1000 | Conductive, Filtered, Sterile | 96pc/rack, 36rack/ctn |
| RTF-50-H-CS | 50 | Filtered, Sterile | 96pc/rack, 36rack/ctn |
| RTF-200-H-CS | 200 | Filtered, Sterile | 96pc/rack, 36rack/ctn |
| RTF-1000-H-CS | 1000 | Filtered, Sterile | 96pc/rack, 36rack/ctn |



Wide Bore Tecan Tips



Beckman Compatible Tips

| Cat No | Max Volume (µL) | Specifications | Packaging |
|---------------|-----------------|-----------------------|---------------------------------|
| RB-50-R-CS | 50 | Non-filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RBF-50-R-CS | 50 | Filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RB-250-R-CS | 250 | Non-filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RBF-50-R-CS | 250 | Filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RB-1000-R-CS | 1000 | Non-filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |
| RBF-1000-R-CS | 1000 | Filtered, Sterile | 96pc/rack, 10rack/box, 5box/ctn |



Hamilton Compatible Tips

| Cat No | Max Volume (µL) | Specifications | Packaging |
|--------------|-----------------|-----------------------------------|-----------------------|
| RH50-H-ES | 50 | Non-filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RHF50-H-ES | 50 | Filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RH50-H-ES | 300 | Non-filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RHF50-H-ES | 300 | Filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RH1000-H-ES | 1000 | Non-filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RHF1000-H-ES | 1000 | Filtered, Sterile, Conductive | 96pc/rack, 36rack/ctn |
| RH50-H-CS | 300 | Non-filtered, Sterile | 96pc/rack, 36rack/ctn |
| RHF50-H-CS | 300 | Non-filtered, Sterile | 96pc/rack, 36rack/ctn |
| RH1000-H-CS | 1000 | Non-filtered, Sterile | 96pc/rack, 36rack/ctn |
| RHF1000-H-CS | 1000 | Non-filtered, Sterile | 96pc/rack, 36rack/ctn |



Enquire for other specifications or packaging



General Consumables

Micro-centrifuge Tubes

Made of transparent polypropylene, the versatile tubes provide a leak proof storage of liquid samples.


- RNase- and DNase- free
- Nonpyrogenic
- Nonsterile tubes are autoclavable
- Graduations
- Low retention specifications available



| Cat No | Volume storage(mL) | Specifications | Max RCF (xg) | Packaging |
|---------------|--------------------|--------------------|--------------|-----------------------|
| MCT-001-060 | 0.6 | Non-Sterile | 12,000 | 1000pc/box, 10box/ctn |
| MCT-001-150 | 1.5 | Non-Sterile | 12,000 | 500pc/box, 10box/ctn |
| MCT-001-150-S | 1.5 | Sterile | 12,000 | 50pc/bag, 5bag/ctn |
| MCT-001-150X | 1.5 | Non-Sterile, Brown | 12,000 | 500pc/box, 10box/ctn |
| MCT-001-175 | 1.7 | Non-Sterile | 12,000 | 500pc/box, 10box/ctn |
| MCT-001-200 | 2.0 | Non-Sterile | 12,000 | 500pc/box, 10box/ctn |
| MCT-001-200S | 2.0 | Sterile | 12,000 | 50pc/bag, 5bag/ctn |
| MCT-001-500 | 5.0 | Non-Sterile | 25,000 | 1050pc/bag, 10bag/ctn |

Enquire for safe lock designed microtubes

Friendly hinge
for finger



Crystal clear for
improved
visibility

High contrast
graduations



Large writing area
for easy
identification



Low Retention Tubes

- RNase- and DNase- free
- Nonsterile
- Autoclavable (121°C, 15min)

| Cat No | Volume storage(mL) | Specifications | Max RCF (xg) | Packaging |
|---------------|--------------------|----------------|--------------|----------------------|
| MCT-001-060LR | 0.6 | Non-Sterile | 12,000 | 500pc/bag, 10bag/ctn |
| MCT-001-150LR | 1.5 | Non-Sterile | 12,000 | 500pc/bag, 10bag/ctn |
| MCT-001-200LR | 2.0 | Sterile | 12,000 | 500pc/bag, 10bag/ctn |





General Consumables

Centrifuge Tubes 15ml | 50mL

Made of transparent polypropylene, the versatile tubes provide a leak proof storage of liquid samples.

- RNase- and DNase- free
- Non-pyrogenic
- Sterile
- Graduations
- Tight sealing screw cap



| Cat No | Volume storage (mL) | Specifications | Max RCF (xg) | Packaging |
|------------|---------------------|----------------|--------------|----------------------|
| CT-002-15A | 15 | Bulk Packed | 15,000 | 25pc/bag, 20bag/ctn |
| CT-012-15A | 15 | Racked | 15,000 | 50pc/rack, 10bag/ctn |
| CT-002-50A | 50 | Bulk Packed | 15,000 | 25pc/bag, 20bag/ctn |
| CT-012-50A | 50 | Racked | 15,000 | 25pc/rack, 00bag/ctn |

Screw Cap Tubes

Made of transparent polypropylene, these colored-cap tubes are suitable for a wide range of temperatures from -80°C to 120°C.

- RNase- and DNase- free
- Non-pyrogenic
- Available with conical or skirted bottom

| Cat No | Volume storage (mL) | Specifications | Packaging |
|---------|---------------------|----------------------|----------------------|
| 801205 | 0.5 | O-ring seal, sterile | 100pc/bag, 10bag/ctn |
| 801207 | 1.5 | O-ring seal, sterile | 100pc/bag, 10bag/ctn |
| 801209 | 2.0 | O-ring seal, sterile | 100pc/bag, 10bag/ctn |
| 801205A | 0.5 | O-ring seal, sterile | 100pc/bag, 10bag/ctn |
| 801207A | 1.5 | O-ring seal, sterile | 100pc/bag, 10bag/ctn |
| 801209A | 2.0 | O-ring seal, sterile | 100pc/bag, 10bag/ctn |

Available Cap Colors







General Consumables


Deep Well Plates

Deep-well plates are used in a wide variety of applications including compound storage, fraction collection, sample mixing and preparation.

- U and V bottom available
- Round and Square mouths
- Autoclavable at 120°C
- Silicone cover mats available



| Cat No | Volume (mL) | Wells | Packaging |
|-----------|-------------|-------|---------------------|
| D-050-96S | 0.5 | 96 | 10pc/bag, 10bag/ctn |
| D-120-96S | 1.2 | 96 | 10pc/bag, 10bag/ctn |
| D-160-96S | 1.6 | 96 | 10pc/bag, 10bag/ctn |
| D-200-96S | 2.0 | 96 | 10pc/bag, 10bag/ctn |
| D-220-96S | 2.2 | 96 | 10pc/bag, 10bag/ctn |
| D-500-48S | 5.0 | 48 | 10pc/bag, 10bag/ctn |
| S-DWP-96 | - | - | 10pc/box, 10box/ctn |
| S-DWP-48 | - | - | 10pc/box, 10box/ctn |






General Consumables

Cryogenic Storage

Manufactured from high-quality virgin polypropylene for storing cells, tissues, microbiological samples, nucleic acid and protein samples in extreme temperatures ranging from -196 to +121°C.


- Tight-sealing and easy opening
- Safe long-term storage down to -196 °C
- Non-pyrogenic
- DNase-and RNase-free
- Sterile



| Cat No | Volume (mL) | Cap Type | Bottom Style | Packaging |
|---------------|-------------|-----------------|--------------|---------------------|
| VC-002-120-EX | 1.2 | External Thread | Round | 50pc/bag, 10bag/ctn |
| CV-002-120-IN | 1.2 | Internal Thread | Round | 50pc/bag, 10bag/ctn |
| VC-002-200-EX | 2.0 | External Thread | Round | 50pc/bag, 10bag/ctn |
| CV-002-200-IN | 2.0 | Internal Thread | Round | 50pc/bag, 10bag/ctn |
| VC-002-500-EX | 5.0 | External Thread | Round | 50pc/bag, 10bag/ctn |
| CV-002-500-IN | 5.0 | Internal Thread | Round | 50pc/bag, 10bag/ctn |

Cryogenic Storage Boxes

| Cat No | Wells | Hinged | Material | Packaging |
|---------|-------|--------|---------------------------|-----------|
| FB-1181 | 81 | No | PP | 100pc/ctn |
| FB-1281 | 81 | Yes | PP | 100pc/ctn |
| FB-1200 | 100 | No | PP | 100pc/ctn |
| FB-1200 | 100 | Yes | PP | 100pc/ctn |
| FB-3281 | 81 | No | Water resistant cardboard | 100pc/ctn |
| FB-3200 | 100 | No | Water resistant cardboard | 100pc/ctn |
| FB-4181 | 81 | Yes | Water resistant cardboard | 100pc/tn |
| FB-4100 | 100 | Yes | Water resistant cardboard | 100pc/ctn |



Hinged boxes available