



Ideal Choice for Digital Dentures

Applications



Denture UV

- ✓ Long-term stability
- ✓ Biocompatible
- ✓ High printing speed
- ✓ MMA free



Denture UV is a biocompatible material developed to produce removable denture bases. This material has excellent process reliability due to its high hardness, and it has low viscosity for easy cleaning and lower material consumption. The unique formulation of the material provides maximum flexural and tensile strength with a stable color tone and odorless characteristics.

Color Pink ●

Specification 1000g Bottle

Certifications FDA

Machine Compatibility



A2D



A2D Ortho



A2D 4K



DS

Properties	Units	Standard	Result
Flexural Strength	MPa	ASTM D790	113.6–128.4
Flexural Modulus	MPa	ASTM D790	2546–2654
Impact Strength (Notched)	J/m	ASTM D648	14.8–19.4
Water Sorption	%	ASTM D570	0.52
Hardness	Shore D	ASTM D2240	87